Allegato "C"

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COVER STAR is the model name for a protective coverall with hood and with elastics at waist, ankles and wrists. **CE-marking** - The coverall corresponds to the European standards of personal protective equipment of category III. The type test and quality assurance certification were issued in 2003 by Sächsisches Textilforschungsinstitut e.V., Annaberger Straße 240, D-09125 Chemnitz (Identification Code 0516). **The European standards for chemical protective clothing** define six types of protection, which are symbolized in one pictogram.

- Typ 1: Gas-tight clothing
- Typ 2: Non-gas-tight clothing
- Typ 3: Liquid-tight clothing
- Typ 4: Spray-tight clothing
- Typ 5: Particle-tight clothing
- Typ 6: Limited splash-tight clothing

AGKI's products specification correspond to the protection types of the European standards. **COVER STAR** coveralls offer protection of Type 5 and Type 6.

The size table combines the body measurements with the standard sizes S-XXXL. Please check your body measurements and select your correct size.

Body measurement (cm)

size	Chest girth	Body height	size	Chest girth	Body height
S	84-92	164-170	XL	108-116	180-188
M	92-100	168-176	XXL	116-124	186-194
L	100-108	174-182	XXXL	124-132	194-202

The symbol of the "Open Book" remembers the wearer to familiarize himself with the "Instructions for Use". **COVER STAR** coveralls are antistatically treated and offer electrostatic protection according to EN 1149-1. In addition **COVER STAR** protection coveralls protection against particulate radioactive contamination according to EN1073-2:2002.

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P=Penetration index (% penetrated) R=Repellency index (%repelled); Remark: For further information concerning the barrier performances please contact

Typical areas for use

COVER STAR coveralls are designed to protect workers against hazardous substances or sensitive products or processes against contamination by man. Depending on the potential risk of the chemical and the circumstances of exposition the **COVER STAR** coveralls can be used for protection against fine particles (Typ5), liquid splashes or liquid spray (Typ6).

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Limitations of use

Please ensure that you have chosen the right garment suitable for your job. **COVER STAR** coveralls offer no penetration resistance against solvents (like n-Heptane and Isopropanol). For advice, please contact your supplier of the garment AGKI®. The user shall be the sole judge for correct combination of full body protective coverall and ancillary equipment (gloves, boots, respiratory protective equipment, etc.) and for how long **COVER STAR** coveralls can be worn with respect to its protective performance, wear comfort and heat stress. AGKI® shall not accept any responsibility whatsoever for improper use of **COVER STAR** coveralls.

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Instructions for use

Marking:

Each Coverall is identified by an inside label. The inside label indicates the type of protection and some further information for use.

1 COVER STAR is the model name for a protective overall with hood and with elastics at waist, ankles and wrists. **2** Manufacturer of the garment. **3** CE-marking - The overall corresponds to the European standards of personal protective equipment of category III. The type test and quality assurance certification were issued in 2003 by Sächsisches Textilforschungsinstitut e.V., Annaberger Straße 240, D-09125 Chemnitz (Identification Code 0516). **4** The European standards for chemical protective clothing define six types of protection, which are symbolized in one pictogram.

- Typ 1: Gas-tight clothing
- Typ 2: Non-gas-tight clothing
- Typ 3: Liquid-tight clothing
- Typ 4: Spray-tight clothing
- Typ 5: Particle-tight clothing
- Typ 6: Limited splash-tight clothing

AGKI's products specification correspond to the protection types of the European standards. COVER STAR coveralls offer protection of Type 5 and Type 6.

5 The size table combines the body measurements with the standard sizes S-XXL. Please check your body measurements and select your correct size.

Body measurement (cm)

size	Chest girth	Body height	size	Chest girth	Body height
S	84-92	164-170	XL	108-116	180-188
M	92-100	168-176	XXL	116-124	186-194
L	100-108	174-182	XXXL	124-132	194-202

6 The symbol of the "Open Book" remembers the wearer to familiarize himself with the "Instructions for Use". **7** COVER STAR coveralls are antistatically treated and offer electrostatic protection according to EN 1149-1. **8** In addition COVER STAR protection coveralls protection against particulate radioactive contamination according to EN1073-2:2002.

9 The five care-pictograms indicate:

- Do not wash. Washing influences the protective performances.
- Inflammable. Keep away from flames or intensive heat.
- Do not machine dry.
- Do not dry clean.
- Do not iron.

Performance profile of COVER STAR coveralls

Physical data	Unit	Result	Test method	Grade
Weight	g/m ²	65		
Abrasion resistance	Cycles	>100; >500	EN530	2
Blocking		No blocking	ISO5978	
Flax cracking resistance	Cycles	>15.000; <40.000	ISO7854	4
Puncture resistance	N	7.63	EN863	1
Mullen burst strength	kPa	71.6	ISO2960	1
Tear resistance	N	Là 62,3 / qu 36,4	ISO9073-4	2
Seam strength	N	108	ISO5082 A2	3
Spraytest Anzug			EN468	erfüllt
Whole-Suit-Test (Staub)			CEN TC162/WG3/N293	erfüllt
Penetration resistance		P R	DIN EN 368	P R
H2SO4 30%		0 98.2		3 3
NaOH 10%		0 98.1		3 3
n-Heptane		0 90.6		3 2
Isopropanol		0 92.8		3 2
Particle penetration (material)	mg	<1	CEN/TC 162/WG3/TG 3 N80	
Protection against radioactive particulates		Nominal protection factor is 15,9	EN1073-2	Class 1

P=Penetration index (% penetrated) R=Repellency index (%repelled);
 Remark: For further information concerning the barrier performances please contact

Typical areas for use

COVER STAR coveralls are designed to protect workers against hazardous substances or sensitive products or processes against contamination by man. Depending on the potential risk of the chemical and the circumstances of exposition the COVER STAR coveralls can be used for protection against fine particles (Typ5), liquid splashes or liquid spray (Typ6).

Preparations before use

Do not use incorrect coveralls. In case of faulty zipper, seams or any other defect, please contact AGKI®. The correct size combined with correct dressing and a closed zipper protected by flap assures the protective performance of the coverall.

Limitations of use

Please ensure that you have chosen the right garment suitable for your job. COVER STAR coveralls offer no penetration resistance against solvents (like n-Heptane and Isopropanol). For advice, please contact your supplier of the garment AGKI®. The user shall be the sole judge for correct combination of full body protective coverall and ancillary equipment (gloves, boots, respiratory protective equipment, etc.) and for how long COVER STAR coveralls can be worn with respect to its protective performance, wear comfort and heat stress. AGKI® shall not accept any responsibility whatsoever for improper use of COVER STAR coveralls.

Storage

COVER STAR coveralls can be stored according to customers storage practices.

Disposal

COVER STAR coveralls can be incinerated or buried in controlled landfill without harm to the environment. Restrictions to disposal result only from contaminants introduced during use.

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